
IKA

designed for scientists

Dispergierwerkzeuge / Dispersion tools / Outils dispersants

S 65 KG-HH-G 65 G / M / F

S 65 KG-HH -D 60 / 80

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Source language: German

Content

	Page
Explication of warning symbols	9
Safety instructions	9
Product information	10
Application	10
Spare parts list	11
Mounting the dispersing tool	11
Cleaning, disinfection and sterilization	11
Maintenance	13
Technical Data	14
Warranty	14

Explication of warning symbols

- DANGER** Indicates an (extreme) hazardous situation, which, if not avoided, will result in death, serious injury.
- WARNING** Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.
- CAUTION** Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
- NOTICE** Indicates practices which, if not avoided, can result in equipment damage.
- CAUTION** Burns hazard!

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

WARNING

 It is only allowed to use these tools on the drive **T65 basic/digital** with the safety switch **SI 400**, which make sure that the drive only starts when the tool is inside the vessel.

WARNING

 Don't touch rotating parts during operation.

CAUTION

 The dispersing tool may heats up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.

CAUTION

 When in operation the dispersing tool can become hot. If the dispersing tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.

CAUTION

 Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There is a risk of splashing and evaporation of liquids; body parts, hair, clothing and jewelry getting caught.

NOTICE

The vessel must be clamped, so that it cannot move during the operation. The **D 60** and **D 80** can destroy the vessel if the dissolver touches the vessel wall. The dissolver tools must have a distance of minimum 100 mm from the vessel wall.

NOTICE

Never run the dispersing tool dry. Without cooling by the medium the sealing will be destroyed.

NOTICE

Rapid temperature change (thermal shock!) can destroy the hard metal sealing surfaces of the dispersing tool **S 65-KG HH**.

- Use the dispersing tool always inside the stirring vessel.
- The work has to be interrupted immediately if you notice unusual noise and/or increasing emission of liquid from the upper, sidewise rinsing drilling. The bearings must be controlled, and replaced if necessary.
- Note the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- Do not use any damaged dispersing tools.
- The dispersing tool is not suitable for continuous operation. Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute. Longer application times bring no improvements, increase however the sample temperature substantially.
- Dispersing causes the medium to heat up.
- Note the min./ max. immersion depth in the “**Technical data**”.
- The function of the dispersing tools depends on the condition of the sharp edges of rotor and stator (integrated in the shaft tube). Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Product information

The dispersing tool consists of a shaft and a generator (stator/rotor). The shaft carries the marking **S 65 KG-HH**, where **KG** stands for the ball bearing and **HH** for a slide ring seal made of hard metal. The seal prevents liquid rising up the shaft and destroying the ball bearing. The top end of the shaft is in the form of a plug-in coupling. The bottom end of the shaft carries a thread for mounting the stator. The shaft – running in ball bearings- carries the rotor.

For the **S 65** range, there are three types of generator, designated **G 65 G**, **G 65 M** and **G 65 F**, where the prefix **G** stands for generator and the suffices **G** for coarse, **M** for medium and **F** for fine. The appropriate generator must be used according to the medium. If the feed material is too coarse the generator can become clogged. The sealing faces then run dry and can be damaged. The generators are suitable for solids components with the following grain sizes:

Generator type	Grain size
G 65 G	40 mm
G 65 M	8 mm
G 65 F	1 mm

It can often be helpful to process a medium with all three generators in sequence, starting with generator **G 65 G**. In such circumstances, to avoid having to replace the entire generator it is worthwhile procuring three shafts, which can be replaced using just a few action steps.

Application

Application: Wet crushing, dispersions, emulsifying (Batch-operation).

Spare parts list

Item	Designation
1	Plastic cap
4	Bearing
9	Coupling
11	O-ring
15	Shaft tube
16	Shaft
18	Rotary seal
19	O-ring
20	O-ring
21	O-ring
30	Stator
40	Rotor


Spare parts diagram see front page.
For orders of spare parts, please specify the type of dispersing tool and the designation of the spare part.

Mounting the dispersing tool

WARNING

The dispersing tool must only be mounted on the dispersion unit when the device is de-energized (e.g. by pulling the plug or pressing the mains switch) and the drive shaft is stationary.

- Open the handwheel bolt (**item 3**) so that the thread does not protrude into the hole.
- Insert the dispersing element (**item 2**) into the drive unit (**item 1**) as far as the stop. The dispersing element is fitted correctly if the recessed edge of the dispersing element is flush with the bottom edge of the drive unit.
- **Attention!** Now secure the dispersing element by screwing the handwheel bolt (**item 3**) tightly.

Cleaning, disinfection and sterilization


Cleaning:

NOTICE

The sterilization of the complete dispersing tool is **not** possible.

The shaft must be cleaned immediately after the unit has been used. This prevents any substance residues adhering to the rotor and stator threads, or bacteria cultures forming in undesired places. Run the dispersion tool in a solvent which will dissolve substance residues but not harm the seals. Because the flow speeds are so high, in such circumstances the rotor and stator are quickly cleaned.


NOTICE

Do not put the dispersion tool in a vessel containing solvent to dissolve substance residues. The seals are mainly dynamic-acting, and quickly start to leak when stationary.

For cleaning purposes disassemble the dispersing tool as described below (items see front page):


Tools for disassembling the dispersing tool:

Note: The tools are included in the delivery of the drive.


Disassemble the rotor and dissolver:

G 65 G:


Insert the shaft wrench between the teeth of the stator. Hold the shaft with the rotor wrench at the coupling section and unscrew the rotor in counterclockwise.

G 65 M/F:


Hold the stator with three pins of the shaft wrench. Hold the shaft with the rotor wrench at the coupling section and unscrew the rotor in counterclockwise.

D 60/80:


Hold the dissolver with a wrench A/F of 17mm (SW17). Hold the shaft with the rotor wrench at the coupling section and unscrew the rotor in counterclockwise.

Note: the user can also hold the shaft with the pin of the rotor wrench and unscrew the rotor with shaft wrench (this procedure will be necessary if the coupling (9) goes loose first.) or SW 17 (for D60/80).

Disassemble the stator:

Hold the shaft tube with hand. Unscrew the stator in counterclockwise (left hand thread) with the hook wrench.

The assembly of the dispersing tool takes place in reverse sequence.

The rotor must not be tightened too strongly, otherwise the thread / rotor teeth will be damaged.

Maintenance

The shaft seals must be regularly checked by dismantling them and examining the contact faces to ensure they are not cracked or pitted. If there is a leak at the seal the suction effect can cause the medium to rise up in the tube around the rotating shaft and may penetrate the drive unit. If the contact faces are cracked or pitted, the rotary seal must be replaced.


Removal and fitting of the rotary seal (see Fig. 1):

To prevent damage to the unit, the slide ring seal must only be changed by a skilled mechanical person.

Change the slide ring seal as follows:

1. Using a A/F 4 hexagon L-wrench to DIN 911, undo the three socket-head threaded pins (3). The lower part of the rotary seal (18b) can then easily be withdrawn.

L-wrench
DIN 911 A/F 4

2. Now the drive shaft can be pushed upwards and out the shaft tube. After this, the upper part of the rotary seal (18a) can be extracted from the shaft tube.
3. When refitting, adjust the bearing pressure of the two sealing surfaces. The correct contact pressure is obtained when the inner part of the lower half of the rotary seal (18c) is pushed in sufficiently far before tightening that it is flush with the outer part of the lower half of the rotary seal (18b).

Spare parts order:

Spare parts diagram see front page.

For orders of spare parts please specify the serial number on the shaft tube, the type of dispersing tool and the designation of the spare part.

Repair:

Please send dispersing tool for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For repair, please request the "Decontamination Certificate" form IKA or use the download printout of it from IKA website: www.ika.com.

If you require servicing, return the instrument in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Technical Data

		S65 KG-HH-G 65 G	S65 KG-HH-G 65 M	S65 KG-HH-G 65 F	S65 KG-HH-D 60	S65 KG-HH-D 80
Working range volumn	l	2 ... 50	2 ... 40	2 ... 30	2 ... 30	2 ... 50
Stator / Rotor (Ø)	mm	65 / 58	65 / 58	65 / 58	- / 60	- / 80
Max. allowable speed	rpm	9500				
Max. circumferential speed	m/s	26.4	28.8	28.8	29.8	39.8
Min. / max. immersion depth	mm	90 ... 450	80 ... 450	80 ... 450	90 ... 450	90 ... 450
Material in contact with medium		FFKM , PTFE, AISI 316L (FFKM and PTFE conform to the FDA, Food and Drug Administration rules)				
pH range		2 ... 13				
Suitable for solvents		yes				
Suitable for abrasive substance		no				
Sterilization methods		all methods				
Working range vacuum	mbar	1				
Temperature range upto	°C	180				
Ultimate fineness dispersions	µm	25 ... 75	20 ... 50	5 ... 20	-	-
Ultimate fieness emulsions	µm	5 ... 25	5 ... 15	1 ... 10	-	-

Subject to technical changes!

Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.
The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited
IKA Works (Thailand) Co. Ltd.