



Viva Dual Access Animal Containment Workstation, Model VDA-4A\_




Viva Universal Animal Containment Workstation, Model VA2-4A\_E



Viva Animal Bedding Disposal Workstation, Model VBD-4A\_

## Animal Containment Workstation

*The Portable Safety Solution for Animal Research Laboratories*



**ESCO**

WORLD CLASS. WORLDWIDE.



## Airflow Sensor

- Monitors real-time airflow for safety
- Alert the user if airflow is insufficient

## Sentinel™ Gold Microprocessor Controller

- Displays all safety information on one screen
- Centered and angled down for easy reach & viewing
- Selectable Quickstart mode for fast operation



## Easy to clean Work Surface and Drain pan

- Two-pieces Stainless Steel Tray, easy to lift.
- Drain hole on both sides to dump animal bedding.



## Easy Work Access

- Large 354 mm (14") Access Opening.
- Accommodates rat and mouse cages.
- Hinged up for easy cleaning.




## Advanced Work Tray Design

- V-shaped Grill to avoid blocking.
- Center Grill to separate Work zone to clean & dirty area.
- Large Tray handle for easy lift.



## Comfortable Leg Room

- 254 mm (10") Leg Room on BOTH sides.
- Reduce fatigue for sitting position.
- Hydraulic Motor to adjust height.



Viva Dual Access Animal Containment Workstation, Model VDA-A\_ Available in 1.2, and 1.5 meter models (4', and 5').

## Accessories and Options

Contact Esco or your Esco Sales Representative for details.

- Electrical Outlets
- Side Shield
- Foldable Side Tray
- Feed Hopper



Side Shield



Feed Hopper





## ELISA Proven Containment


- Provides >99% Allergen Containment.
- Ensures User's Safety.



VIVA®

## ULPA Filter

- 10x Filtration efficiency of HEPA filter
- Creates ISO Class 3 work zone instead of industry-standard ISO Class 5




## Quiet Operation

- The quietest Dual-Access Animal Workstation in the world, at 53 dbA in open field condition
- Comfortable for the operator and animals




## ISOCIDE™ Powder Coat

- Silver-ion impregnated powder coat
- Inhibit microbial growth to improve safety



## Dual Energy Efficient DC ECM Motor

- Powered by latest generation DC ECM motor, that is more efficient than legacy ECM and VFD motors
- 70% Energy savings compared to AC motor
- Stable airflow, despite building voltage fluctuations & filter loading



| Standards Compliance | Air Quality   | Filtration  | Electrical Safety  |
|----------------------|---|---|--|
|                      | ISO 14644.1, Class 4, Worldwide<br>JIS B9920, Class 4, Japan<br>JIS B55295, Class 4, Japan<br>US Fed Std 209E, Class 10 USA | EN-1822 (H14), Europe<br>IEST-RP-CC001.3, USA<br>IEST-RP-CC007, USA<br>IEST-RP-CC034.1, USA | UL-61010A-1, USA<br>CSA22.2, No.1010-192, Canada<br>EN61010-1, Europe<br>IEC61010-1, International |

**ESCO**

WORLD CLASS. WORLDWIDE.



**ESCO**

Escp UA2  
Air-flow: OK  
I: 90 fpm  
Socket: ON

15:34  
Sash: OK  
D: 70 fpm

**Sentinel<sup>™</sup>GOLD**

MICROPROCESSOR CONTROL SYSTEM

MENU

VSET  
DIAGNOSTICS


## Sentinel<sup>™</sup> Gold Microprocessor Controller

- Displays all safety information on one screen
- Centered and angled down for easy reach & viewing
- Selectable Quickstart mode for fast operation



## Single-Piece Wall

- Large radius for easy cleaning
- Side-mounted electrical outlets and staggered service fixtures, for easy reach



## Single-Piece Work Tray

- Recessed to contain spillage
- Curved grill to prevent blockage



## Raised Arm Rest

- Helps prevent grille blocking
- Comfortable working posture



## Angled Drain Pan

- Helps prevent grille blocking
- Does not harbor contaminants

## Airflow Sensor

- Monitors real-time airflow for safety
- Alert the user if airflow is insufficient



Viva Universal Animal Containment  
Workstation, Model VA2-\_A-E

Available in 1.2, and 1.8 meter models (4' and 6').

**VIVA**


Lab Animal Research Products • Universal Animal Containment Workstation

**EN12469  
COMPLIANT**




## Energy Efficient DC ECM Motor

- Powered by latest generation DC ECM motor, that is more efficient than legacy ECM and VFD motors
- 70% Energy savings compared to AC motor
- Stable airflow, despite building voltage fluctuations & filter loading




## ULPA Filter

- 10x Filtration efficiency of HEPA filter
- Creates ISO Class 3 work zone instead of industry-standard ISO Class 5



## ISOCIDE™ Powder Coat

- Silver-ion impregnated powder coat
- Inhibit microbial growth to improve safety



## ELISA Proven Containment

- Provides >99% Allergen Containment.
- Ensures User's Safety.



| Standards Compliance | Air Quality  | Filtration  | Electrical Safety  |
|----------------------|--|---|--|
|                      | ISO 14644.1, Class 3, Worldwide<br>JIS B9920, Class 3, Japan<br>JIS B55295, Class 3, Japan<br>US Fed Std 209E, Class 1 USA | EN-1822 (H14), Europe<br>IEST-RP-CC001.3, USA<br>IEST-RP-CC007, USA<br>IEST-RP-CC034.1, USA | UL-61010A-1, USA<br>CSA22.2, No.1010-192, Canada<br>EN61010-1, Europe<br>IEC61010-1, International |



### Sentinel™ Silver Microprocessor Controller

- Displays all safety information on one screen
- Centered and angled down for easy reach & viewing
- Selectable Quickstart mode for fast operation



### Bang Bars

- Increase efficiency of bedding disposal operations.



### Integrated Waste Chute

- Dispose refuse bag safely within the work zone



#### Operator and Environmental Protection

The VIVA Bedding Disposal Workstation provides operator and environmental protection from animal allergen.

#### Exclusive hydraulic height-adjustable stand

Allows the work surface height to be adjusted to user preference therefore minimizing strain during repetitive operations.

### Airflow Sensor

- Monitors real-time airflow for safety
- Alert the user if airflow is insufficient



Viva Animal Bedding Disposal Workstation Model VBD-4A\_

Available in 1.2 meter model (4') only.






### Carbon Filter

- Nanocarb™ activated carbon filter removes odors



### ULPA Filter

- 10x Filtration efficiency of HEPA filter
- Yields 10x cleaner lab air from allergen than industry standard HEPA filter




### ELISA Proven Containment

- Provides >99% Allergen Containment.
- Ensures User's Safety.




### ISOCIDE™ Powder Coat

- Silver-ion impregnated powder coat
- Inhibit microbial growth to improve safety



| Standards Compliance | Filtration  | Electrical Safety |
|----------------------|---|-------------------|
|                      | EN-1822 (H14), Europe<br>IEST-RP-CC001.3, USA<br>IEST-RP-CC007, USA<br>IEST-RP-CC034.1, USA | UL61010-1, USA    |






### VDA Cabinet Airflow System

- The VDA Dual Access Workstation employs recirculating airflow configuration for better filtration efficiency.
- The blower system pulls ambient intake air through the front grilles, creating inflow that provides operator protection from allergen inside the work zone. An activated carbon pre-filter removes odors.
- Air flows through the common plenum on top of the cabinet. A portion of it goes up through ULPA filter as exhaust to create inflow. The remaining portion goes down through ULPA supply filter and bathes the work zone in clean air with a non-turbulent downflow.

- ULPA-filtered air
- Unfiltered / Potentially contaminated air
- Room air / Inflow air

- The combination of vertical laminar downflow and inflow creates an air curtain to protect the operator from contaminants released from the work surface.




### VA2 Cabinet Airflow System

- Ambient air pulled through the perforations towards the work zone front prevents contamination of the work surface and work product. The inflow does not mix with the clean air within the cabinet work zone. Inflow air travels through a return path toward the common air plenum (blower plenum) at the top of the cabinet.
- Approximately 40% of the air in the common plenum is exhausted through the ULPA filter to the room. The remaining 60% of the air is passed through the downflow ULPA filter and into the work area as a vertical laminar flow air stream bathing the work surface in clean air.

- ULPA-filtered air
- Unfiltered / Potentially contaminated air
- Room air / Inflow air

- The uniform, non-turbulent air stream protects against cross-contamination within and throughout the work area.
- Near the work surface, the ULPA-filtered downflow air stream splits with a portion moving toward the front air grille, and the remainder moving to the rear air grille. A small portion of the downflow enters the side capture zones at a higher velocity (small blue arrows).
- A combination of inflow and downflow air streams form an air barrier that prevents contaminated room air from entering the work zone, and prevents work surface emissions from escaping the work zone.



### VBD Cabinet Airflow System

- Carbon Filter
- Blower
- Exhaust ULPA Filter
- Pre-Filter

- Room air is drawn in across the front of the cabinet with an average velocity of 0.35 m/s (70 fpm).
- Air is drawn up through the cabinet's work zone and forced through the ULPA filter (>99.999% typical efficiency for 0.1 to 0.3 micron sized particles).

- ULPA-filtered air
- Unfiltered / Potentially contaminated air
- Room air / Inflow air

- The full work zone ceiling extraction system ensures airflow uniformity throughout the cabinet's main chamber.
- The ULPA filtered air then returns to the laboratory stripped of all airborne contaminants and odor.






## General Specifications, VIVA® Dual Access Animal Containment Workstation, Model VDA

| Model                                      |                                   | VDA-4A_  | VDA-5A_  |
|--|-----------------------------------|--|--|
| External Dimensions (W x D x H)            |                                   | 1340 x 762 x 1961 mm (52.8" x 30.0" x 77.2") min height<br>1340 x 762 x 2245 mm (52.8" x 30.0" x 88.4") max height | 1645 x 762 x 1961 mm (64.7" x 30.0" x 77.2") min height<br>1645 x 762 x 2245 mm (64.7" x 30.0" x 88.4") max height |
| Internal Work Area (W x D x H)             |                                   | 1100 x 465 x 564 mm (43.3" x 18.3" x 22.2")  | 1405 x 465 x 564 mm (55.3" x 18.3" x 22.2")  |
| Downflow Velocity                          |                                   | 0.24 m/s (47 fpm)  |  |
| Pre-Filter                                 |                                   | Disposable and non-washable polyester fibres with 85% arrestence / EU3 rated                                       |  |
| ULPA Filter Typical Efficiency             |                                   | >99.999% for particle size between 0.1 to 0.3 microns, per IEST-RP-CC001.3   |  |
| Sound Emission per EN 12469*               |                                   | 53 dBA   | 54 dBA   |
| Fluorescent Lamp Intensity at Zero Ambient |                                   | 1725 lux (160 foot candles)  | 1525 lux (142 foot candles)  |
| Construction, Main Body                    |                                   | 1.5 mm (0.06") 16 gauge EG Steel with Isocide™ Oven-Baked Epoxy-Polyester Powder Coated Finish                     |  |
| Shipping Dimensions, Maximum (W x D x H)   |                                   | 1720 x 820 x 2240 mm (67.7" x 32.2" x 88.1")   | 2025 x 820 x 2240 mm (79.7" x 32.2" x 88.1")   |
| Shipping Weight                            |                                   | 342 Kg (754 lbs)   | 432 Kg (952 lbs)   |
| Shipping Volume, Maximum                   |                                   | 3.16 m³ (111.6 cu.ft.)   | 3.72 m³ (131.4 cu.ft.)   |
| Electrical Rating                          | VDA-_A8                           | 220-240 VAC, 50 / 60 Hz, 1Ø  |  |
|  | VDA-_A9                           | 110-130 VAC, 50 / 60 Hz, 1Ø  |  |
| Power Consumption                          | VDA-_A8                           | 190 W  | 230 W  |
|  | VDA-_A9                           | 210 W  | 250 W  |
| Accessories                                | Foldable Side Tray (SS Shelf Kit) | VDA-001 5170257  |  |
|  | Side Shield                       | VDA-004 5170562  | VDA-005 5170563  |
|  | Feed Hopper                       | VDA-006 5170594  |  |

\* Noise as measured in open field / anechoic chamber.



**ESCO**


WORLD CLASS. WORLDWIDE.

## General Specifications, VIVA® Universal Animal Containment Workstation, Model VA2

Note to customer: Insert electrical voltage number into last model number digit\_when ordering.

| Model  |               | VA2-4A_ -E  |  | VA2-6A_ -E                                      |                                      |
|--|---------------|---|--|---|--------------------------------------|
| Nominal Size   |               | 1.2 meter (4')  |  | 1.8 meter (6')                                  |                                      |
| External Dimensions (W x D x H)                                  |               | 1423 x 815 x 1510 mm<br>(56" x 32.1" x 59.4")   |  | 2030 x 815 x 1510 mm<br>(79.9" x 32.1" x 59.4") |                                      |
| Maximum External Dimensions with Support Stand (W x D x H)       |               | 1585 x 852 x 2235 mm<br>(62.4" x 33.5" x 88.0")   |  | 2193 x 852 x 2235 mm<br>(86.3" x 33.5" x 88.0") |                                      |
| Internal Work Area (W x D x H)                                   |               | 1270 x 623 x 680 mm<br>(50.0" x 24.5" x 26.7")  |  | 1870 x 620 x 680 mm<br>(73.6" x 24.4" x 26.7")  |                                      |
| Average Airflow Velocity   | Inflow        | 0.45 m/s (90 fpm)   |  |   |                                      |
|  | Downflow      | 0.35 m/s (70 fpm)   |  |   |                                      |
| Airflow Volume   | Inflow        | 625 m³ / h (368 cfm)  |  | 921 m³ / h (542 cfm)                            |                                      |
|  | Downflow, 60% | 959 m³ / h (547 cfm)  |  | 1414 m³ / h (832 cfm)                           |                                      |
|  | Exhaust, 40%  | 625 m³ / h (368 cfm)  |  | 921 m³ / h (542 cfm)                            |                                      |
| ULPA Filter Typical Efficiency                                   |               | >99.999% for particle size between 0.1 to 0.3 microns per IEST-RP-CC001.3                   |  |   |                                      |
| Sound Emission*  | NSF / ANSI 49 | 63 dBA  |  | 64 dBA  |                                      |
|  | EN 12469      | 60 dBA  |  | 61 dBA  |                                      |
| Fluorescent Lamp Intensity                                       |               | > 1400 lux (> 130 foot candles)   |  | > 1230 lux (> 114 foot candles)                 |                                      |
| Cabinet Construction   |               | 1.5 mm (16 gauge) electrogalvanized steel with Isocide white oven-baked epoxy power coating |  |   |                                      |
| Net Weight Cabinet including stand                               |               | 406 Kg (895 lbs)  |  | 528 Kg (1164 lbs)                               |                                      |
| Shipping Weight Cabinet including stand                          |               | 456 Kg (1005 lbs)   |  | 570 Kg (1257 lbs)                               |                                      |
| Shipping Dimensions, Maximum (W x D x H) Cabinet excluding stand |               | 1550 x 950 x 1900 mm (61.0" x 37.4" x 74.8")  |  | 2150 x 950 x 1900 mm (84.6" x 37.4" x 74.8")    |                                      |
| Shipping Volume, excluding stand                                 |               | 2.80 m³ (99 cu.ft.)   |  | 3.88 m³ (137 cu.ft.)                            |                                      |
| Electrical*  | Model         | Voltage   |  | Model   | Voltage                              |
|  | VA2-4A1-E     | 220-240 VAC, 50/60 Hz, 1Ph, 5.5 amps  |  | VA2-6A1-E                                       | 220-240V, AC, 50/60 Hz, 1Ph, 6 amps  |
|  | VA2-4A2-E     | 110-120 VAC, 50/60 Hz, 1Ph, 11 amps   |  | VA2-6A2-E                                       | 110-120V, AC, 50/60 Hz, 1Ph, 12 amps |

\* Noise as measured in open field / anechoic chamber.



1. Airflow Sensor
2. Exhaust ULPA / H14 Filter
3. Energy-efficient DC ECM Blower
4. Downflow ULPA / H14 Filter
5. UV light Retrofit Kit Provision
6. IV bar Retrofit Kit Provision
7. Pre-filter

8. RS 232 Port, zero volt relay contacts for exhaust and alarm system
9. Electrical Panel
10. Fluorescent Lamps
11. Plugged Service Fixture provisions (2 on each side)
12. Electrical Outlet Retrofit Kit Provision
13. Stainless Steel Single-piece Work Tray

14. Stainless Steel Arm Rest
15. Drain Valve Retrofit Kit Provision
16. Sentinel™ Gold Microprocessor Control System
17. Safety Glass Sliding Sash Window
18. Single-piece Stainless Steel Back Wall and Side Walls
19. Removable Side Panel for plumbing access




## General Specifications, VIVA® Bedding Disposal Workstation, Model VBD-4A\_

| Nominal Size                             |                              | 1.2 meter (4')   |                        |                        |
|--|------------------------------|--|------------------------|------------------------|
| External Dimensions (W x D x H)          |                              | 1247 x 760 x 1966 mm (49.1" x 30.0" x 77.4") minimum height<br>1247 x 760 x 2271 mm (49.1" x 30.0" x 89.4") maximum height |                        |                        |
| Internal Work Area (W x D x H)           |                              | 1040 x 680 x 594 mm (40.9" x 26.8" x 23.4")  |                        |                        |
| Work Surface Height                      |                              | 920 mm ~ 1225 mm (36.2" ~ 48.2")   |                        |                        |
| Front Opening                            |                              | 400 mm (15.7")   |                        |                        |
| Inflow Velocity                          |                              | 0.35 m/s (70 fpm) at initial setpoint  |                        |                        |
| Pre-Filter                               |                              | Disposable, non-washable polyester fiber, 85% arrestance, EU3 rated  |                        |                        |
| ULPA Filter Typical Efficiency           |                              | >99.999% at 0.1 to 0.3 microns as per IEST-RP-CC001.3 USA  |                        |                        |
| Sound Emission* Per EN 12469             |                              | 58 dBA   |                        |                        |
| Fluorescent Lamps                        |                              | > 1,300 lux (> 121 foot candles)   |                        |                        |
| Workstation Construction                 | Main Body                    | 1.2 mm (0.05") 18 gauge electro-galvanized steel with Isocide™ white oven-baked epoxy-polyester powder-coating             |                        |                        |
|  | Work Top                     | 1.2 mm (0.05") 18 gauge stainless steel, type 304, with 4B finish  |                        |                        |
|  | Inner Liner                  | 0.9 mm (0.035") 20 gauge stainless steel, type 304, with 4B finish   |                        |                        |
| Net Weight                               |                              | 233 Kg (514 lbs)   |                        |                        |
| Shipping Weight                          |                              | 294 Kg (648 lbs)   |                        |                        |
| Shipping Dimensions, Maximum (W x D x H) |                              | 2150 x 1840 x 1230 mm (84.6" x 72.4" x 48.4")  |                        |                        |
| Shipping Volume, Maximum                 |                              | 4.87 m³ (172 cu.ft.)   |                        |                        |
| Electrical**                             | Model                        | VBD-4A1  | VBD-4A2                | VBD-4A3                |
|  | Voltages                     | 220-240 VAC, 50 Hz, 1Φ   | 110-120 VAC, 60 Hz, 1Φ | 220-240 VAC, 60 Hz, 1Φ |
|  | Cabinet Full Load Amps (FLA) | 3 A  | 6.5 A                  | 3 A                    |
|  | Optional Outlets FLA         | 5 A  | 5 A                    | 5 A                    |
|  | Cabinet Nominal Power        | 309 W  | 268 W                  | 309 W                  |
|  | Cabinet BTU                  | 1054   | 914                    | 1054                   |

11


\* Noise as measured in open field / anechoic chamber.



**ESCO**

WORLD CLASS. WORLDWIDE.

## ESCO GLOBAL NETWORK



ART Equipment  
Biological Safety Cabinets  
CO<sub>2</sub> Incubators  
Compounding Pharmacy Equipment  
Containment / Pharma Products  
Ductless Fume Hoods  
Freeze Dryer  
Lab Animal Research Products  
Laboratory Fume Hoods  
Laboratory Ovens and Incubators  
Laminar Flow Clean Benches  
PCR Cabinets  
PCR Thermal Cyclers  
Powder Weighing Balance Enclosures  
Ultra-low Freezers

*The Esco Group of Companies is a global life sciences tools provider with sales in over 100 countries. The group is active in lab equipment, pharma equipment and medical devices. Manufacturing facilities are located in Asia and Europe. R&D is conducted worldwide spanning the US, Europe and Asia. Sales, service and marketing subsidiaries are located in 12 major markets including the US, UK, Singapore, Japan, China and India. Regional distribution centers are located in the US, UK, and Singapore.*

Life Science • Chemical Research • Assisted Reproductive  
Technology (ART) • Pharmaceutical Equipment • General Equipment

# ESCO

WORLD CLASS. WORLDWIDE.

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777  
Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com  
www.escoglobal.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA  
Toll-Free USA and Canada 1-877-479-3726 • Tel 215-441-9661 • Fax 484-698-7757  
eti.sales@escoglobal.com • www.escolifesciences.us

**Esco Global Offices:** Bahrain | Bangladesh | China | India | Indonesia | Italy | Japan | Malaysia  
Philippines | Russia | Singapore | South Africa | South Korea | Thailand | United Kingdom | USA | Vietnam

**dijkstra**  
**vereenigde**

LABORATORIUM-APPARATUUR.NL  
T. +31 (0)320-266 171 | laboratorium@dijkstra.net  
pascallaan 9 - 8218 NJ Lelystad | www.dijkstra.net

escoglobal.com  
LIVE SUPPORT

ISOCIDE™

