

Traditional hydrolysis apparatus



Manual acid hydrolysis for up to six samples. Followed by classic fat extraction according to Weibull-Stoldt with **traditional Soxhlet** units or automatic **SOXTHERM** extraction systems.

Complete unit, ready for operation: with special supports, glass condensers, ring shaped rest, silicone rubber tubings, beakers (1000 ml), filter insert with 6 funnels, 100 folded filters, soft bottle with narrow mouth (1000 ml), 250 g boiling stones and power cord.

Available models

- HY 16/4 Hydrolysis unit, 4-place, 230 V 50-60 Hz, Order no. 10-0019
- HY 16/6 Hydrolysis unit, 6-place, 230 V 50-60 Hz, Order no. 10-0020

For information on our fully automatic hydrolysing system, please refer to **HYDROTHERM**.



Postbus 2151 Pascallaan 9
8203 AD Lelystad 8218 NJ Lelystad
Tel: 0320-266171 Fax: 0320-257354
email: laboratorium@dijkstra.net
www.dijkstra.net