

# GALAXY 230 COLONY COUNTER



## Features

- ▶ Innovation PC software makes your data traceable
- ▶ Pressure sensor design enables the counter being operated by any marking pen
- ▶ Wide range size of petri dish holder
- ▶ Background plate changeable (white, black)
- ▶ Unique BACK button allow you to deduct data when over counting
- ▶ Non-glare illumination
- ▶ Adjustable pressure sensitivity and acoustical counting signal

## Description

*galaxy 230* is designed for reliable and traceable counting of bacterial colonies growing on nutrient agar in petri dishes with a wide range of diameter from 60 to 150mm.

Innovation PC software empower you recording, tracing, or reporting data of sample No., date, user's name, row count, dilution factor, total count etc.

Pressure sensor design enables the counter being operated by any marking pen

Ring-shaped illumination lamp with electronic ballast results in glare-free light without blinding or noise and makes operator working comfortable without eyestrain

Unique BACK button allow you to deduct data when over counting Contrast-plate with white and black side is supplied for using dark or clear agar.

## Specification

Readout:	4-Digit (0-9999) bright red 0.6" LED display
Magnification:	2~3x
Petri dish size:	90-150mm diameter
Zero reset:	Yes
Count back:	Yes
Illumination:	Glare-free ring-shaped lamp
Pressure sensitivity:	Adjustable
Acoustical counting signal:	Adjustable
Out put:	RS232
Dimension:	W313mm x H360mm x D346 mm
Weight:	Net: 5.2 KG Shipping: 6.4 KG



Innovation PC software empower you recording, tracing, or reporting data of sample No., date, user's name, row count, dilution factor, total count etc.

## Ordering Information

175200-11	galaxy 230 colony counter, Power AC110V, 60Hz, complete with marking pen (Black) 1 piece
175200-22	Same as 175200-11, but AC220V, 50Hz
175300-11	galaxy 230 colony counter, Power AC110V, 60Hz, complete with marking pen (Black) 1 piece, PC software
175300-22	Same as 175300-11 but AC220V, 50Hz
175200-48	STAEDTLER 4 color (black, red, blue, green) marking pen set, 0.4mm, permanent ink
175300-12	PC software with connecting cable for Galaxy 230 colony counter



Postbus 2151                      Pascallaan 9  
8203 AD Lelystad              8218 NJ Lelystad  
Tel: 0320-266171              Fax: 0320-257354  
email: [laboratorium@dijkstra.net](mailto:laboratorium@dijkstra.net)  
[www.dijkstra.net](http://www.dijkstra.net)

**NEW**

### Dragon

Electronic Instant Ignition Burner



Electronic Instant Ignition burners are also available for sterilizing of Manifolds, Platinum loops etc.