

SONOSHAKE

NEW

The new generation: Ultrasonic bath SONOREX DIGITEC DT 1028 F combined with shaking device SA 1028



Ultrasonic bath DT 1028 F + shaking device SA 1028 = SONOSHAKE

Code No: 3257

SONOREX DIGITEC DT 1028 F

- ⇒ especially for sonication of samples in laboratory flasks and for fast degassing
- ⇒ flat ultrasonic bath with a stainless steel tank, inner dimensions 500 × 300 × 65 mm (l × w × d)
- ⇒ housing made of stainless steel, protected against spray, 535 × 325 × 205 mm (l × w × h)
- ⇒ outlet: ball valve G ½
- ⇒ ultrasonic peak output 1280 W
- ⇒ HF output 320 W
- ⇒ working frequency 35 kHz
- ⇒ Current consumption 1,4 A

⇒ required floor space of SONOSHAKE approx. 850 × 360 mm (l × w)

The SONOSHAKE can be used either separately as an ultrasonic bath for sample preparation, a shaking device SA 1028 only, or both combined!

Code No.: 3249

SA 1028 attachable to existing SONOREX DIGITEC DT 1028 F

SA 1028

- ⇒ analogue setting of time and shaking frequency
- ⇒ reciprocating motion: settings in 4 steps possible of up to 200/min
- ⇒ constant amplitude of 20 mm independently of loading
- ⇒ rack easy to remove – easy mounting of the laboratory clamps EK 10 – 100 (separately to order)
- ⇒ floor space approx. 410 × 280 mm (l × w)
- ⇒ mounting of 36 × 10-ml-flasks or 36 × 25-ml-flasks or 18 × 50-ml-flasks or 12 × 100-ml-flasks

